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Conventional electrodes are normally single junction. As depicted by the ?figure above, these electrodes have only a single junction between the internal reference wire and the external solution. Under adverse conditions, such as high pressure, high temperature, highly acidic or alkaline solutions, the positive flow of the electrolyte through the junction is often reversed resulting in the ingress of sample solution into the reference compartment. If this is left unchecked, the reference electrode can become contaminated, leading to complete electrode failure. Another potential problem with single junction electrodes is the clogging of the junction due to silver chloride (AgCl) precipitation. Silver can be easily precipitate in samples that contain Tris buffer or heavy metals. When the electrolyte solution makes contact with the sample, some AgCl will precipitate on the external face of the junction. The result is drift readings obtained from the sensor.

Hannaâ€™s double junction system, as the name implies, has two junctions, only one of which is in contact with the sample as shown in the ?figure. Under adverse conditions, the same tendency of sample ingress is evident. However, as the reference electrode system is separated physically from the intermediate electrolyte area, the contamination of the electrode is minimized. The likelihood?of clogging of the junction is also reduced with a double junction electrode since the outer reference cell uses a ?fill solution that is ?silver-free?. Since there is no silver present, there is no precipitate that can form to clog the junction.

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El cuerpo del HI1332B estÃ¡ compuesto de resina de polieterimida (PEI). El cuerpo de PEI es adecuado para una amplia gama de aplicaciones y destaca en las mediciones de terreno debido a su durabilidad. El blindaje alrededor de la punta estÃ©rica de vidrio minimiza la rotura debido a golpes o caÃ±adas accidentales del electrodo. El plÃ¡stico PEI es un plÃ¡stico de alta calidad que es quÃ³micamente resistente a muchos productos quÃ³micos agresivos.